




KING CAPSULE

Expert in Empty Hard Capsules



HPMC * PULLULAN * GELATIN



Overview



KING CAPS is a leading brand for empty hard capsules in China. The factory full name Zhejiang Pujiang Enerkang Capsule Co., Ltd. KING CAPS focuses on production of HPMC capsules, pullulan capsules and gelatin capsules. The long-term commitment allows us to be able to provide professional solution based on customized demands.

KING CAPS was founded by a group of experienced experts in R&D, production and QC of capsules. In order to maintain quality level, we adopted advanced automatic capsule production lines. The workshop is designed and constructed in strict accordance with GMP standard. Product & quality system are certified by USDA-Organic/GMP/ FDA/ISO/HALAL/KOSHER.

Most commonly used capsule sizes are available: 000#, 00#, 0#, 1#, 2#, 3#, 4#.

Meanwhile, KING CAPS provides a variety of colors for choice. We can also produce according to pantone color or customer's sample. Our team will be pleased to work with you on your unique solution! you on your unique solution!



Our Certifications



Quality Control

Good Manufacturing Practices (GMP's) are the guidelines that provide the procedures, processes and documentation to insure the product manufacturing process consistently meets the expected specifications:

- Identity
- Potency
- Composition
- Product Purity

KING CAPS ensures quality begins with the incoming receipt of raw materials and packaging components, through WIP (work in process) to the shipment of all finished products.

Raw material quality control:

In order to guarantee safety and stability, all the raw materials for our capsule is inspected rigidly according to Enterprise quality standards which is made on the bases of State pharmacopoeia.

Quality Control in Manufacturing:

- Sampling inspection is conducted during the producing process to keep all the parameters in right normal scope to strengthen quality stability.
- Glue dissolution: The temperature, color and viscosity are monitored under strict control.
- Capsule production: appearance, dimension and moisture are checked by quality controller. Environment temperature and moisture level in factory room, as well as glue viscosity and temperature are monitored by the operator.
- Screening Process: The sorted capsule by sorting machine is inspected by quality controller further to ensure all the capsule for sale meet quality standards. The sorting parameter is adjusted and controlled by workers in charge of sorting machine running.

Finished Product Management

We have the latest equipment to test inbound raw material for efficacy, pureness, and potency. We also emplo the latest test equipment to ensure your product is free from contaminants and matches all label claims when it leaves our facility.

Each batch record and formulation are stored in our secure database which is backed up on a daily basis. Records are maintained on-site as well as off-site to guarantee security in the event of a disaster.

Physical & chemical indexes of finished capsule would be checked before packaging.
Only the qualified capsule would be packed for clients.

Tracking of Product:

According to GMP regulation, all the capsule are packed with clear batch number and the quality records are made for tracking if necessary.

Quality Management System:

Product & quality system are certified by USDA-Organic/GMP/FDA/ISO/HALAL/KOSHER.



HPMC Empty Capsules

Advantages

HPMC (Hydroxypropyl methylcellulose) Empty Capsules are two-piece capsules made from cellulosic raw materials that satisfy vegetarian, religion, cultural and policy needs. These vegetable capsules are an attractive, all-natural dosage form that retain all the advantages - easy to swallow, effectively mask taste and odor, and allow product visibility. HPMC capsules are also starch-free, gluten-free and preservative-free, and meet the strict dietary needs of customers that choose a vegetarian lifestyle. HPMC capsules have been developed for both Health Supplements products and Nutraceutical.



Features

- Non-animal origin, allergen free, gluten-free and preservative-free.
- Free of amino acids, no risk of cross-linking reactions & Maillard reaction.
- Low moisture content and high stability, ideal capsule for high hygroscopic contents.
- Easy to be disintegrated in short time to release the stuffing inside.

Worldwide Acceptance

KING CAPS HPMC empty capsules meet the requirements of Islam, Judaism and vegetarians, also meets the needs of a wide variety of health-conscious consumers.

HPMC has been recorded in pharmacopoeia worldwide, such as American USP, European EP, Japanese. It has been widely adopted in the production and packaging of medicine, cosmetics and healthcare supplements.



Pullulan Empty Capsules



Advantages

KING CAPS's pullulan empty capsules, a vegetarian capsule made from naturally fermented pullulan for a more extensive and distinctive appeal to the booming and discerning healthy lifestyle market worldwide.

Pullulan capsules are the premium capsules designed for the needs of the growing and very discerning healthy lifestyle market.

Pullulan capsules provide the premier balance of performance and purity for healthy lifestyle consumers who are looking for naturally healthy attributes in the products they buy. Made from tapioca which is naturally fermented into pullulan, a starch-free vegetarian capsule, Pullulan capsules speak to the needs of the most discerning consumers.

Features

- No toxic chemicals involved. Completely natural and environmentally-friendly
- Starch-free, gluten-free, preservative-free, vegetable origin and non-GMO
- Excellent oxygen barrier properties
- Strong visual appeal, brilliantly clear and visually appealing to customers
- Highly compatible with encapsulation machines

Sizes & Specifications

KING CAPS offers various colors and printing styles for Pullulan capsules. We can also manufacture with customized dimensions and colors. Our professional team will be pleased to work with you on your unique solution!

Size		Length	Lock length	Pre-Lock length	Everage Weight	Volume	Reference volume based on powder density, in mg				Wall thickness	Outer diameter
		±0.20mm	±0.3mm	±0.5mm	mg	ml	0.6g/ml	0.8g/ml	1.0g/ml	1.2g/ml	±0.02mm	mm
000#	cap	12.6	26.2	28.0	200±8	1.4	840	1120	1400	1680	0.115~0.135	10.00~10.10
	body	22.2										
00#	cap	11.7	23.3	25.3	140±8	0.95	570	760	950	1140	0.110~0.130	8.50~8.60
	body	20.0										
0#	cap	11.0	21.6	23.6	108±6	0.68	400	533	667	800	0.105~0.125	7.61~7.71
	body	18.6										
1#	cap	9.8	19.3	21.3	83±5	0.5	300	400	500	600	0.100~0.120	6.90~7.00
	body	16.7										
2#	cap	9	17.7	19.7	68±4	0.37	220	293	367	440	0.080~0.100	6.32~6.40
	body	15.5										
3#	cap	8.1	15.8	17.8	55±4	0.3	170	227	283	340	0.080~0.100	6.32~6.40
	body	13.6										
4#	cap	7.2	14.1	16.1	44±3	0.21	120	160	200	240	0.080~0.100	6.32~6.40
	body	12.3										



Zhejiang Pujiang Enerkang Capsule Co., Ltd.
Add: No.588, Fenghe Road, Pujiang, Zhejiang, China
Tel: +86 (0532) 87680032 +86 18562696939
Web: www.kingcapsule.com
E-mail: info@kingcapsule.com

